

Specialty Mattresses

For Hospital and Patient Care bed




- **Non-allergic**
- **Durable, safe materials**
- **Easy-to-clean, hygienic and anti-bacterial**
- **Soft feel, firm support**
- **Versatile range to meet user demands**


The range of Merivaara specialty mattresses has been developed to meet the unique demands of a variety of patient care environments. They have been designed for use on Merivaara beds, with the best qualities of each product being utilized to ensure optimal patient comfort. Various quality options are available to meet different application needs. These mattresses use only the highest quality materials that have been thoroughly tested under laboratory and actual patient care conditions.

Specialty mattresses


For hospital and patient care beds




Medi-Form 10000 0948 200 × 75 × 13 cm
10000 0949 200 × 85 × 13 cm



Medi-Ultra 10000 0931 200 × 75 × 12 cm
10000 0933 200 × 85 × 12 cm
12800 9933 205 × 80 × 12 cm
12300 9913 166 × 80 × 12 cm



Medi-Flex-Safe 10000 0941 200 × 75 × 13 cm
10000 0943 200 × 85 × 13 cm



Medi-Flex-Tricot 10000 0937 200 × 75 × 13 cm
10000 0938 200 × 85 × 13 cm



Medi-Tricot 10000 0935 200 × 75 × 10 cm
10000 0936 200 × 85 × 10 cm

Medi-Form

Extremely soft mattress with profiled surface for institutional use, mattress forms according to temperature and patient body

- prevents pressure sores and numbing
- best suited for long-term use, even in higher patient weights
- flexible, conforms easily to adjustable beds
- top layer: made of visco-elastic Formax foam
base layer: Relax-55 foam
- cover: easy-to-clean, bi-elastic, hygienic anti-bacterial fabric
- fire classification: SL1

Medi-Ultra

High quality foam mattress intended for extended care needs

- soft top, sturdy, hygienic mattress with dual-sided profile
- the soft top **reduces pressure** and prevents patient numbing
- a durable mattress suited for a range of patient weights
- flexible, conforms well to adjustable beds
- made of high elastic Relax-55 foam
- cover: easy-to-clean, bi-elastic, hygienic anti-bacterial fabric
- fire classification: SL1

Medi-Flex-Safe

Very flexible form-cut hygienic mattress

- best suited for short-term use for persons of medium weight
- form cutting allows the mattress to bend in both directions
- made of E-33 PS high resiliency technical foam
- cover: easy-to-maintain, bi-elastic, hygienic anti-bacterial fabric
- fire classification: SL1

Medi-Flex-Tricot

Standard form-cut mattress for institutional use; recommended for use on 3- or 4-sectional beds requiring frequent base position adjustments

- best suited for short-term use for patients of medium weight
- form cutting allows the mattress to bend in both directions
- made of E-35 resilient foam
- tricot cover
- fire classification: SL2

Medi-Tricot

Firm standard mattress for institutional use

- best suited for short-term use for patients of medium weight
- made of E-35 resilient foam
- tricot cover
- fire classification: SL2



Merivaara Corp.

Puustellintie 2, FI-15150 LAHTI, FINLAND
Tel. +358 3 3394 611 • Fax +358 3 3394 6144
merivaara@merivaara.fi

www.merivaara.com

Washing instructions clearly marked on mattress covers.