

Merivaara Today



Merivaara Corp.



Merivaara has been a specialist manufacturer of high-quality medical equipment for more than 100 years, turning out its first operating table in 1910. The company manufactures all core equipment for today's modern hospitals. It covers a wide spectrum of hospital furniture, including operating tables, operating and examination lights, beds and bed stands, delivery beds, stretchers, and examination tables. Products developed and manufactured by Merivaara are distributed worldwide to more than 100 countries.

All products are CE labelled and Merivaara's quality system complies with the ISO 9001:2000 and ISO 13485:2003 standards. The company's environmental management system is ISO 14001:2004 certified.

Merivaara Range of the Operating Tables



Promerix



A heavy duty operating table with excellent performance and features and modular table top - easy to set up for all procedures

- Heavy duty mobile operating table - max. safe working load 325 kg / sliding version 275 kg
- Modular construction - Easy detachable table sections
- Available control systems: cable hand control, IR- hand control, override panel, foot control
- Pre-programmed positions + customer programmable positions
- Fulfills X-ray and C-arm requirements

Practico



Electric mobile operating table for general surgery with modern stylish design and easy detachable table top sections

- Modular operating table for general and day surgery
- With separate attachments suitable for i.e. gyn/ urology, eye/ENT and even orthopaedics incl. shoulder arthroscopy procedures - Beach Chair
- Max. lifting capacity 380 kg
- Safe working load 225 kg
- Optional 300 mm horizontal sliding
- Battery-operated with integrated charger/ transformer, optional main power use
- Back-up control for emergency situations
- Wide range of accessories



Rapido

Needs of the entire day surgery process is covered with one product as Rapido operates as a transport trolley, operating table and recovery trolley

- Open design allows easy access to the patient
- Easy to manoeuvre
- Wide table top with a soft mattress gives extra comfort during recovery
- Full X-ray capability in each table section
- Combines economy and versatility as the table can accommodate multiple applications
- Also available is a model for ENT surgery
- Wide range of accessories



Emergo

Emergo provides advanced features and an ergonomic design. This modular construction stretcher can be used for patient transport, first aid, trauma cases, recovery and minor procedures

- Excellent ergonomic features, e.g. subframe and drop-down pushbars
- 2- or 3- sectional X-ray imaging platform and back section X-ray rack base for trauma patients
- Safe working load is 225 kg
- Reliable and safe in treatment situations requiring high speed
- Modular construction
- Fifth wheel for light and quick transport



Optima

Upgraded model for different delivery positions, enhanced ergonomics for midwife and mother

- A simultaneously sliding and rising movement in back section reduces pressure in abdominal area
- Numerous comfortable resting, labour and delivery positions (e.g. squatting and sitting positions)
- Ergonomic working positions for midwives - leg section gets down and the bed height can be adjusted up to 970 mm
- Safe working load 230 kg, for leg section 150 kg
- Safety for mother and infant: wide mattressbase and foldable siderails
- Several colour options

Merilux / Meriled

Merivaara offers both halogen and led lights for hospitals. Meriled light series for surgical procedures demanding highest optical quality and Merilux light family with triple IR-filtering optical systems for examination and surgery



Merilux

- Merilux X1, X3 and X5 can be combined to Duo or Trio
- Merilux Vision has X3 and/or X5 with camera and/or screen option
- Exceptional depth of illumination
- Outstanding laminar flow characteristics
- Easy and quick bulb changes without tools
- Easy to clean, and sterilize
- Shadowless and steady light output

Meriled

- Cool light
- Long effective led life 30 000 - 50 000 hours
- Low power consumption 150 W
- Wireless controls
- Endlight (green)
- No colored shadows
- Uniform shape of light field
- Light combinations: Solo, Duo 360, Trio
- Light intensity 160.000 lux (1 m)
- High intensity illumination 160 klux
- Easy to manoeuvre



Carena

Carena, the next-generation hospital bed is equally suited to emergency room and long-term care, intensive care or nursing homes

- High patient weight capacity 275 kg
- Unique lower height 295 mm and upper height 890 mm
- Modular construction, innovative structure, new design
- Unparalleled patient ergonomics
- State-of-the art ergonomics for the nurse
- Long service life-cycle and ease of the maintenance



Merivaara Range of the Hospital Beds



Carena

Adatto

Futura Plus

Child Cot

Child Cot

Futura Plus Jr

Futura Plus

Adatto

Carena



MERIVAARA

Merivaara Corp.

Puustellintie 2, FI-15150 LAHTI, FINLAND

Tel. +358 3 3394 611, fax +358 3 3394 6144

merivaara@merivaara.fi

www.merivaara.com