

Ada / Nova

Bedside cabinets



Ada

The most easy-to-use reading and dinner table adjuster mechanism on the market:

- can be adjusted with one hand
- stepless height adjustment
- handedness can be easily and quickly changed without any tools



The twin castors make the cabinet easy to move. Castor locks anchor the cabinet firmly in place, a drawer lock is standard equipment. The cabinet is also available without the reading/dinner table (Figure a). With a wide range of accessories, the cabinet can be modified to meet specific needs.

The Ada Bedside Cabinet frame is made of epoxy-coated steel, top and side panels are light grey, birch or beech melamine. There are three colour options for the drawer fronts. The plastic drawers are lightweight and easy to clean.

Technical specifications and accessories

Dimensions

- width: 475 mm (dinner table width: 545 mm)
- height: 900 mm
- depth: 600–1130 mm (with dinner table in front)
- castor Ø: 75 mm (individually-locking)

Side panels

- grey melamine.....10000 9186
- birch melamine10000 9187
- beech melamine.....10000 9188

Drawers

- white fronts10000 9210
- light blue fronts10000 9212
- blue fronts.....10000 9218

Dinner tables

- reading and dinner table.....10000 9190

Accessories

- table top cover ABS plastic (b)10000 9262
- magazine basket (f).....10000 9266
- slipper rack (b,d)10000 9268
- patient intercom unit holder10000 0629
- bottle holder for lower drawer (e).....10000 9270
- tray for small items, narrow (e).....10000 9280
- tray for small items, wide (f).....10000 9285
- bottle or waste bag rack (g).....10000 9288

warranty 3 years



Nova



2670



2660

Nova Bedside Cabinet 2670 and 2660 have easy-to-clean surfaces. Twin-castors ensure easy and quiet movability.



2660

Nova Bedside Cabinet including reading/dinner table. The balanced construction makes the bedside cabinet safe when using the reading/dinner table. Adjustments can be conveniently made regardless of the bed's position.



2781

Separate holder for patient intercom unit can be attached with finger screw. Colours according to the drawers.



2680

Nova Reading/Dinner Table 2680. The Nova collection includes a separate stand-mounted dinner/reading table.

Nova – a stylish bedside cabinets are functional and efficient to hospital and nursing home care facilities. Modern materials and manufacturing technology make the bedside cabinets easy to handle and keep clean. The sturdy, adjustable reading/dinner table can be accessed from either side of the cabinet. The easy gliding castors ensure that the bedside cabinet moves easily and quietly. Rounded edges and reliable functions ensure safe use.

Technical specifications and accessories

- Frame of anodized aluminium profile and composite material
- Drawers made of melamine or wooden laminate (max. load 50 kg)
- Twin castors Ø 50 mm/75 mm, separate castor locks
- The height of reading/dinner table steplessly adjustable with gas springs
- Reading/dinner table made of composite material (max. load 10 kg)

Dimensions

- Width: 2660: 660 mm (reading/dinner table lowered)
2670: 480mm
2680: 790 mm/620 mm (reading/dinner table lowered)
- Depth 2660: 880 mm (75mm Ø castors)
2670: 780 mm (50mm Ø castors)
- Height 2660: 880 mm (75mm Ø castors)
2670: 780 mm (50mm Ø castors)
- Reading/dinner tray (model 2660 and 2680)
520 × 360 mm (max. load 10 kg)
- Adjustment range of the reading/dinner table (model 2660 and 2680) 795–1075 mm with Ø 75 mm castors

Accessories and options

- table top cover, ABS-plastic 120002762
- holder for patient intercom unit:
 - stained rose-wood laminate 12000 2794
 - stained beech laminate 12000 2792
 - birch laminate 12000 2790
 - light grey melamine 12000 2781
- cabinet colour options:
 - stained rose-wood laminate
 - stained beech laminate
 - birch laminate
 - light grey melamine

warranty 3 years



A stylish bedside cabinet series



MERIVAARA

Merivaara Corp.

Puustellintie 2, FI-15150 LAHTI, FINLAND

Tel. +358 3 3394 611, fax +358 3 3394 6144

merivaara@merivaara.fi

www.merivaara.com