

Carro de mantenimiento 285

Maintenace trolley for operating rooms.

- Stainless steel frame
- Fixed surfaces in stainless steel, raised edges prevent things from falling out
- Push bar
- Castors (ø 125 mm) in corrosion resistant material, two locking castors
- Upper level height from floor: 800 mm



Datos Técnicos

Ancho	560 mm
Longitud	800 mm
Altura	900 mm
Tamaño de las ruedas	125
Código de producto	100000285

Contacta con nosotras

Merivaara Corp.

Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com