

# Head rest 18101

Head rest 18101 is designed for Rapido operating table for general surgery, but it fits to all Merivaara operating tables with Beach Chair back rest when you need to perform other surgical procedures. Highly recommended for Smarter Practico Beach Chair. Requires no adapter for Promerix and Practico operating tables with Beach Chair back rest. For Standard and Kidney bridge back rests the fixing adapter 60146 is needed. Use product code 145181 when ordering the head rest with the adapter.

## Order information:

145181 Head rest 18101 with adapter 60146 100  
018 101 Headrest 100 060 146 Adapter for  
Promerix and Practico



## Technical Specifications

Width	408 mm
Height	300 mm
Thickness	80 mm
Max. patient weight	325 kg
Safe working load (SWL)	14 kg
Trendelenburg	-45 °
Reverse Trendelenburg	+25 °
Weight	4.1 kg
Product Code	145181, 100018101

## **Contact us**

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)