

ARM REST 125

Padded arm board 125 includes a fixing clamp and straps 110.
Adjustable height, gas spring assisted angle adjustment. Rotates
180 degrees. Suitable for EUR-rails and USA-rails.

Order no.

110125 Arm rest 125 incl. clamp for EUR, US or UK rails



Width	640 mm
Height	205 mm
Thickness	60 mm
Max. patient weight	325 kg
Safe working load (SWL)	19.5 kg
Trendelenburg	-36 °
Reverse Trendelenburg	+26 °
Weight	6.2 kg
Product Code	110125