

# THIGH HOLDER I2012

Rotating, complete with antistatic pad. Vertical bar 550mm,  
horizontal bar 184mm.

Includes universal clamp I0200 for EUR, UK and US rails.



Order no.

100 012 012 Thigh holder I2012 (1pc), includes universal clamp I0200 for EUR, UK and US rails

<b>Length</b>	100 mm
<b>Height</b>	550 mm
<b>Max. patient weight</b>	325 kg
<b>Safe working load (SWL)</b>	40 kg
<b>Weight</b>	1.8 kg
<b>Product Code</b>	I00012012