

Merivaara screen protector specification

Technical specification

Code	Description	Width (mm)	Height (mm)	Corner radius (mm)	Material	Compatible monitor
A43817300	Display protector, rounded, 654.8x422, R17.5 mm	654,8	422	17.5	Transparant Medical 5H	Foreseeson FM-A2701D, FM-A2701DS
A43817400	Display protector, view area, 698.4x392.9, R0 mm	698,4	392,9	0	Transparant Medical 5H	Foreseeson FS-L3202D
A43817500	Display protector, rounded, 648x400, R5 mm	648	400	5	Transparant Medical 5H	Foreseeson FS-L2702DT, FM-B2702D, FS-L2702D
A43817700	Display protector, rounded, 748x428, R5 mm	748	428	5	Transparant Medical 5H	Foreseeson FS-P3102DNB
A43817800	Display protector, rounded, 748x453, R8 mm	748	453	8	Transparant Medical 5H	Foreseeson FM-E3203D
A43817900	Display protector, rounded, 593x412, R19 mm	593	412	19	Transparant Medical 5H	Barco MDSC-2324
A43818100	Display protector, rounded, 638x417, R19 mm	638	417	19	Transparant Medical 5H	Barco MDSC-2326
A43818200	Display protector, rounded, 767x462, R14 mm	767	462	14	Transparant Medical 5H	Barco MDSC-2232, MDSC-8231
A43818300	Display protector, view area, 648x356, R0 mm	648	356	0	Transparant Medical 5H	Barco MDSC-8427

Installation

Installation instructions: T406000