


Disposable sterile light handle cover


Disposable sterile light handle cover for Q-Flow surgical light and examination light Merilux X1

- Sterile lamp handles protection for single use
- Protection size fits various types of handles
- Lamp handle with protection 150 mm allows access the light switch on the handle, if applicable
- The perforated plastic collar acts both as plain protection against light and while also flexible holder
- The installed plastic cap creates a barrier between the operator and handle
- Sterile packed in peel pack, Ethylene oxide sterilization

A43460601 Set of 200 units

Merivaara Corp. 04.11.2025





Technical Specifications

Product Code	A43460600

Merivaara Corp. 04.11.2025



Contacta con nosotras

Merivaara Corp.

Tarmontie 2-4 15860 HOLLOLA FINLAND

+358 3 3394 611

merivaara@merivaara.com

Merivaara Corp. 04.11.2025