

## Tube mount for the handle of Q-Flow Fluent

- Adjustable
- Material: steel
- Includes 1 pc handle
- Compatible with detachable sterilizable handle 485492 for Q-Flow Fluent
- To be used with Vista monitor holder



### Technical Specifications

Product Code	A43932300

### Contacta con nosotras

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)