

# TRIO

## TRIO suspension solution for surgical lights and monitors

Modern operating rooms today are equipped with more and more equipment designed to provide surgical staff both the clinical images as well as critical video sources during the procedure. For the staff, functionality and ergonomics of suspension solutions that carry these vital monitors, is of very high importance.

Suspension solutions of Q-Flow product family has gone through a complete update, now being able to offer many new solutions that will enhance the ergonomics of any operating room.

[Read how Merivaara's R&D was challenged to think outside the box and what then happened.](#)

### EXAMPLE OF AVAILABLE ARM OPTION:

- **TRIO** 110 mm ceiling tube



## Technical Specifications

Product Code	520230

## Contacta con nosotras

**Merivaara Corp.**

Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)