

## Hand fixation device 10393

- Weinberger-type hand fixation device
- For hand and finger traction
- Can be used only for hand operating table 10390
- Translucent and sterilizable.
- Size: 180 mm x 179 mm
- Weight: 0,3 kg



### Technical Specifications

Weight	0,3 kg
Product code	10393, 100010393

### Contacta con nosotras

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)