

## Traction device 11443

Traction device 11443 for intramedullary nailing of the humerus. To be used with Weinberger hand fixation device 10393.



### Technical Specifications

Width	676 mm
Length	932 mm
Weight	3 kg
Adjustment range	0-200 mm
Product code	11443, 100011443

## Contacta con nosotras

**Merivaara Corp.**

Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)