

# DORO Teflon Head Rest System

- Autoclavable and machine washable system which conforms to WHO Infection Control Guidelines for cleaning and sterilization
- Suitable also for high risk infection patients!
- 3-pin rigid cranial fixation
- For supine, prone, lateral or sitting position (for sitting position Crossbar Adaptor is required)
- Excellent stability and flexibility

## Includes the items below:

103 003 009      *Skull Clamp Teflon*  
103 002 009      *Doro Swivel Adaptor Teflon*  
103 001 009      *Doro Adjustable Base Unit Teflo.*



## Technical Specifications

Product Code	100304000
--------------	-----------

## Contacta con nosotras

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)

