

Head rest double articulation, narrow HR1120

- versatile adjustment offers optimum patient positioning for surgery
- angular adjustments with gas springs
- adjustment angle +37°/-58° and height adjustment 135 mm
- requires adapter 60146
- especially suitable for ENT
- Can be used without adapter for ENT Upper-Back Section 60325 and EYE ENT145325



TECHNICAL SPECIFICATIONS

| | |
|--------------|--------|
| | |
| Product Code | HR1120 |

Contacta con nosotras

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com