


Oxygen bottle holder 18189

For holding 5-10 L oxygen bottle on the Rapido operating tables's base cover.

- For bottle diameter ca. 140 mm
- Max bottle length with valves 550 mm
- To be fitted by the customer, mounting instructions included
- Nylon fabric straps with zinc plated locking buckles
- Suitable with Rapido



Technical Specifications

Вес	2.25 kg
код продукта	100018189

Contacta con nosotras

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com