

Multi lateral support, wide SS1050

- multi-adjustable lateral support for optimum positioning
- padding (50 mm) designed to prevent pressure and nerve injuries
- antimicrobial tightening levers
- height adjustment range 200 mm
- pad size 200x200 mm
- Recommended clamps: CL2010/ CL2110



TECHNICAL SPECIFICATIONS

Product code	SS1050

Contacta con nosotras

Merivaara Corp.

Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com