

TrenGuard HYBRID -potilastuki Trendelenburg-asentoon

TrenGuard on edistyksellinen ratkaisu potilaan asetteluun MIS- ja robottiaivusteisessa kirurgiassa.

Hyödyt:

- Ehkäisee potilaan liukumista Trendelenburg-asennossa erityisen niskakuoppapehmusteen avulla.
- Ei tarvetta olkatuille. Pehmeät lateraaliset stabilointityynyntä kontrolloivat vartalon painopisteen muutoksia.
- Ei tarvetta vartalohihnoille helpottaen potilaan ventilaatiota.
- Tehokas ja helppokäyttöinen. Sopii kaikkiin Merivaaran leikkauspöytiin (myydään vain leikkauspöydän mukana).
- Vähemmän jätettä uudelleen käytettävien lateraalisten stabilointityynyntä ansiosta.
- Alhaisemmat materiaalikustannukset.
- Järjestelmä sopii jopa 250 kg painoisille potilaille.



450 Hybrid starter package Standard

Including:

- Patient Support Frame
- Clamps for EUR rails
- Single-use soft head stabilizing pillow and Speed Bump bolster
- Lateral Stabilizing Pillows
- Product code A43824200

Reusable package

Including:

- One reusable soft lateral stabilizing pillow control body mass shift during transition to Trendelenburg
- Order No. 55251

Reusable package for bariatric patients

Including:

- One reusable soft lateral stabilizing pillow to control body mass shift during transition to Trendelenburg
- Order No. 56251

Single-use wedge package

For bariatric patients and those without discernible cervical notch.

Including:

- 12 single-use integrated wedge and head stabilizing pillows
- Order No. 56401

Single-use package

Including:

- 12 single-use Speed Bump bolsters which prevent patient sliding
- 12 single-use soft head stabilizing pillows which control side-to-side rocking of the head
- Order No. 55201

Single-use package for bariatric patients

Including:

- 12 single-use Speed Bump bolsters which prevent patient sliding
- 12 single-use soft head stabilizing pillows which control side-to-side rocking of the head
- Order No. 56201

Ota yhteyttä

Merivaara Corp.

Tarmontie 2-4
15860 HOLLOLA
SUOMI

+358 3 3394 611

merivaara@merivaara.com