

K-Tek 250 SAAPASMALLINEN JALKATUKI

The Kyra 250 is a junior stirrup designed for patients aged 6-12 years old but also effective in smaller/lighter adults weighing up to 250 lbs/113 kg.

- More working space with wider boot rotational range
- Lightweight and Strong with Kyra's Core-Grip™ technology

Order information

100012814	K-Tek 250 stirrup pair
100012808	K-Tek™ stirrup cart
100012811	K-Tek™ stirrup boot pad kit
100 011 043	Rail clamp, EUR and UK
100 011 043U	Rail clamp, US



Tekniset Tiedot

Potilaan enimm. paino	113 kg
Paino	4.6 kg
Toimintasäde	+ -25
Vertikaalinen säätöalue	+85 - -55 °
Abduktioalueen laajuus	+25 - -9 °
Tuotekoodi	100012814

Ota yhteyttä

Merivaara Corp.

Tarmontie 2-4
15860 HOLLOLA
SUOMI

+358 3 3394 611

merivaara@merivaara.com