

K-Tek 350 SAAPASMALLINEN JALKATUKI

The K-Tek Comfort Stirrups for patient positioning for a wide range of gynecologic, urologic and laparoscopic procedures. The Stirrups increase surgeons' operating room performance by permitting simple adjustment of lithotomy and abduction intra-operatively. The K-tek Comfort 350 stirrup is designed for patients with a BMI of less than 29 with a large flexible lateral fin and can support patients up to 159 kg weight capacity.

- Designed to relieve pressure with soft padding that encapsulates the foot.
- K-TEK™ silicone closure straps making cleaning easy.
- Flexible finned boots accommodates various leg sizes.
- Superior access with lithotomy range of motion from -55° to +85° aids in relieving clinician shoulder/arm strain.



Order information

100012806	K-Tek Comfort 350 stirrups
100012808	K-Tek™ stirrup cart
100012811	K-Tek™ Comfort™ 350 stirrup additional boot pad kit (one piece)
CL1210	Pikalukko, lattavarsi
100 011 043U	Rail clamp, US



Tekniset Tiedot

Potilaan enimm. paino	159 kg
Lithotomy Range	+85 - -55 °
Abduction range	+35 - -9 °
Product Code	100012806

Ota yhteyttä

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
SUOMI

+358 3 3394 611

merivaara@merivaara.com