

# K-Tek 500 SAAPASMALLINEN JALKATUKI

The K-Tek Comfort Stirrups for patient positioning for a wide range of gynecologic, urologic and laparoscopic procedures. The Stirrups increase surgeons' operating room performance by permitting simple adjustment of lithotomy and abduction intra-operatively. The K-tek Comfort 500 stirrup is designed for patients with a BMI of greater than 29 with a large flexible lateral fin and can support patients up to 227 kg weight capacity.

- Designed to relieve pressure with soft padding that encapsulates the foot.
- K-TEK™ silicone closure straps making cleaning easy.
- Flexible finned boots accomodates various leg sizes.
- Superior access with lithotomy range of motion from -55° to +85° aids in relieving clinican shoulder/arm strain.



## Order information

100012812	K-Tek 500 stirrup pair
100012808	K-Tek™ stirrup cart
100012811	K-Tek™ stirrup boot pad kit
100 011 043	Rail clamp, EUR and UK
100 011 043U	Rail clamp, US

## Tekniset Tiedot

Potilaan enimm. paino	227 kg
Paino	6.6 kg

Toimintasäde	+ -25
Lithotomy Range	+85 - -55 °
Abduction range	+25 - -9 °
Product Code	100012812

## Ota yhteyttä

**Merivaara Corp.**

Tarmontie 2-4  
15860 HOLLOLA  
SUOMI

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)