

## K-Tek 800 stirrup pair 12813

The K-Tek Comfort Stirrups for patient positioning for a wide range of gynecologic, urologic and laparoscopic procedures. The Stirrups increase surgeons' operating room performance by permitting simple adjustment of lithotomy and abduction intra-operatively. The K-tek Comfort 800 stirrup is designed for patients with a BMI of greater than 29 with a large flexible lateral fin and can support patients up to 363 kg weight capacity. The boot has wider boot design to facilitate a larger mass.

- Designed to relieve pressure with soft padding that encapsulates the foot.
- K-TEK™ silicone closure straps making cleaning easy.
- Flexible finned boots accommodates various leg sizes.
- Superior access with lithotomy range of motion from -55° to +85° aids in relieving clinician shoulder/arm strain.



### Order information

100012813	K-Tek 800 stirrup pair
100012808	K-Tek™ stirrup cart
100012811	K-Tek™ stirrup boot pad kit
CL1210	Pikalukko, lattavarsi
100 011 043U	Rail clamp, US

### Technical Specifications

Max. patient weight	363 kg
---------------------	--------

Weight	6.6 kg
Movement angle	+25
Lithotomy Range	+85 - -55 °
Abduction range	+25 - -9 °
Product Code	100012813

## Ota yhteyttä

**Merivaara Corp.**

Tarmontie 2-4  
15860 HOLLOLA  
SUOMI

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)