K-Tek 800 stirrup pair 12813

The K-Tek Comfort Stirrups for patient positioning for a wide range of gynecologic, urologic and laparoscopic procedures. The Stirrups increase surgeons' operating room performance by permitting simple adjustment of lithotomy and abduction intra-operatively. The K-tek Comfort 800 stirrup is designed for patients with a BMI of greater than 29 with a large flexible lateral fin and can support patients up to 363 kg weight capacity. The boot has wider boot design to facilitate a larger mass.

- Designed to relieve pressure with soft padding that encapsulates the foot.
- K-TEK™ silicone closure straps making cleaning easy.
- Flexible finned boots accomodates various leg sizes.
- Superior access with lithotomy range of motion from -55° to +85° aids in relieving clinican shoulder/arm strain.



Order information

100012813 K-Tek 800

stirrup pair

100012808 K-Tek™

stirrup cart

100012811 K-Tek™

stirrup boot

pad kit

CL1210 Pikalukko,

lattavarsi

100 011 043U Rail clamp,

US

Technical Specifications

Max. patient weight	363 kg

Merivaara Corp. 19.11.2025



Weight	6.6 kg
Movement angle	+-25
Lithotomy Range	+8555 °
Abduction range	+259 °
Product Code	100012813

Ota yhteyttä

Merivaara Corp. Tarmontie 2-4 15860 HOLLOLA SUOMI

+358 3 3394 611

merivaara@merivaara.com

Merivaara Corp. 19.11.2025