

K-Tek 800 stirrup pair 12813

The K-Tek Comfort Stirrups for patient positioning for a wide range of gynecologic, urologic and laparoscopic procedures. The Stirrups increase surgeons' operating room performance by permitting simple adjustment of lithotomy and abduction intra-operatively. The K-tek Comfort 800 stirrup is designed for patients with a BMI of greater than 29 with a large flexible lateral fin and can support patients up to 363 kg weight capacity. The boot has wider boot design to facilitate a larger mass.

- Designed to relieve pressure with soft padding that encapsulates the foot.
- K-TEK™ silicone closure straps making cleaning easy.
- Flexible finned boots accomodates various leg sizes.
- Superior access with lithotomy range of motion from -55° to +85° aids in relieving clinican shoulder/arm strain.



Order information

| | |
|--------------|-----------------------------|
| 100012813 | K-Tek 800 stirrup pair |
| 100012808 | K-Tek™ stirrup cart |
| 100012811 | K-Tek™ stirrup boot pad kit |
| CL1210 | Pikalukko, lattavarsi |
| 100 011 043U | Rail clamp, US |

Technical Specifications

| | |
|---------------------|--------|
| Max. patient weight | 363 kg |
|---------------------|--------|

| | |
|-----------------|-------------|
| Weight | 6.6 kg |
| Movement angle | +25 |
| Lithotomy Range | +85 - -55 ° |
| Abduction range | +25 - -9 ° |
| Product Code | 100012813 |

Ota yhteyttä

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
SUOMI

+358 3 3394 611

merivaara@merivaara.com