

Kevyt jalkaosa moottorijalkaosiin 60204

Light leg section improves work ergonomics of the operating room staff. It is very light and easy to handle. Light weight leg section attaches to the operating table with Fluent locking mechanism. Light Leg Section for motorized legs, 100060204, is suitable to be used with Promerix operating table and New, smarter Practico with motorized leg joints. Includes 80 mm thick, antistatic VEF (viscoelastic foam) mattress.



Technical Specifications

Длина	630 mm
ширина (столешница)	540 mm
Высота	76 mm
Безопасная рабочая нагрузка	40 kg
Вес	4.9 kg
код продукта	100060204

Ota yhteyttä

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
SUOMI

+358 3 3394 611

merivaara@merivaara.com