

Kyra 400 SAAPASMALLINEN JALKATUKI

The K-tek 400 stirrup is designed for bariatric patients.

- Non-finned boot
- Can be used with robotic surgery
- Easy to clean durable silicone straps
- Includes more padding to reducing pressure
- Increased range of motion



Order information

100 012 815	K-Tek 400 stirrup pair
100012808	K-Tek™ stirrup cart
100012811	K-Tek™ stirrup boot pad kit
100 011 043	Rail clamp, EUR and UK
100 011 043U	Rail clamp, US

Tekniset Tiedot

Макс. вес пациента	181 kg
код продукта	100012815

Ota yhteyttä

Merivaara Corp.

Tarmontie 2-4
15860 HOLLOLA
SUOMI

+358 3 3394 611

merivaara@merivaara.com