

Light non-motorized leg section 60203

Light leg section improves work ergonomics of the operating room staff. It is very light and easy to handle. Light weight leg section attaches to the operating table with Fluent locking mechanism. Light leg section for non-motorized legs, 100060203, is suitable for Smarter Practico operating table with non-motorized leg joints, Practico electro-mechanical and Practico Manual operating tables. To be used with Promerix operating table, 25- or 40- section is needed first to be attached to the seat section. Includes 80 mm thick, antistatic VEF (viscoelastic foam) mattress.



Tekniset tiedot

Pituus	630 mm
Leveys (makuutaso)	540 mm
Korkeus	76 mm
Turvallinen työkuorma	40 kg
Paino	4.9 kg
Tuotekoodi	100060203

Ota yhteyttä

Merivaara Corp.

Tarmontie 2-4
15860 HOLLOLA
SUOMI

+358 3 3394 611

merivaara@merivaara.com