

K-Tek Arm Boards

Integrated clamp snap locks onto rail at any position. Designed to prevent fluid ingress. Radiolucent, ultradurable strong material. Eliminates gear grinding and pinch hazards. One-hand adjustable. Order no. 100012818 K-TEK2005 Arm Board, EU-rail 100012819 K-TEK2000 Arm Board, US-rail 100012820 K-TEK2006 Arm Board, UK-rail



Technical Specifications

Max. patient weight	227 kg
Weight	1.8 kg
Product code	100012818, 100012819, 100012820

Ota yhteyttä

Merivaara Corp.

Tarmontie 2-4
15860 HOLLOLA
SUOMI

+358 3 3394 611

merivaara@merivaara.com