

K-Tek2340 Secondary Shoulder Support Pair

K-Tek2340 Secondary Shoulder Supports are designed to provide non-contact support for extreme Trendelenburg positioning. Flexible swivel adjustment for optimal positioning. Pads included. Sold in pairs. Requires K-TEK4800 Round Post Rail Clamps.



Technical Specifications

Max. patient weight	272 kg
Product Code	100012823

Ota yhteyttä

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
SUOMI

+358 3 3394 611

merivaara@merivaara.com