

## K-Tek2340 Secondary Shoulder Support Pair

K-Tek2340 Secondary Shoulder Supports are designed to provide non-contact support for extreme Trendelenburg positioning. Flexible swivel adjustment for optimal positioning. Pads included. Sold in pairs. Requires K-TEK4800 Round Post Rail Clamps.



### Technical Specifications

Max. patient weight	272 kg
Product Code	100012823

## Ota yhteyttä

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
SUOMI

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)