

# Infusion rod 444

Mobile IV-pole with pushing handle.

- Height adjustment: a. 1510-2210
- IV-pole material: Stainless steel
- 5-wheel
- Wheel diameter  $\varnothing$  80 mm
- Stand diameter a. 680 mm
- 4-hooks
- Color of plastic parts: Grey RAL 7035
- Material of handle: Stainless steel



## Technical Specifications

Product code	100000444
--------------	-----------

### Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)

