

Instrument Trolley 609

Instrument trolley for operating rooms.

- Chromed steel frame
- 2 removable stainless steel trays
- Castors \varnothing 75 mm for easy moving
- Trolleys fit under an 80 cm high working level



Technical Specifications

Width	400 mm
Length	600 mm
Height	750 mm
Castors size	75
Product code	100000609

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com