

3D monitors

Merivaara offers a wide range of medically certified monitors suitable for visualizing the content during medical procedures. We offer also a range of compatible accessories and even customized mounting solutions to meet your demands.



Available models

SUPPLIER CODE

EX3220-3D

EX3141-3D

FM-E3204DGC

FM-A5505DGC

PRODUCT CODE

A45000735

A43818700

A43861300

A43861400

Contact us

Merivaara Corp.

Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com