

4K/UHD displays for monitor arms

Merivaara offers a wide range of medically certified monitors suitable for visualizing the content during medical procedures. We offer also a range of compatible accessories and even customized mounting solutions to meet your demands.



Available models

SUPPLIER CODE

[FM-B2702DG \(12G-SDI\)](#)
[FM-E3203D](#)
[FM-E3203DG](#)
[MDSC-8527 SSTP](#)
[MDSC-8532 SSTP](#)

PRODUCT CODE

A43861100
A43824500
A43816900
TBD
TBD

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com