


## 4K/UHD displays for monitor arms

Merivaara offers a wide range of medically certified monitors suitable for visualizing the content during medical procedures. We offer also a range of compatible accessories and even customized mounting solutions to meet your demands.



### Available models

#### **SUPPLIER CODE**

[FM-B2702DG \(12G-SDI\)](#)  
[FM-E3203D](#)  
[FM-E3203DG](#)  
[MDSC-8527 SSTP](#)  
[MDSC-8532 SSTP](#)

#### **PRODUCT CODE**

A43861100  
A43824500  
A43816900  
TBD  
TBD

## Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)