

Other Full HD monitors

Merivaara offers a wide range of medically certified monitors suitable for visualizing the content during medical procedures. We offer also a range of compatible accessories and even customized mounting solutions to meet your demands.



Available models

SUPPLIER CODE

[FS-E2101D](#)

FS-P2404D

FS-P2604D

PRODUCT CODE

A43818400

A45000469

A43966000

Contact us

Merivaara Corp.

Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com