


## Q-Flow compatible Full HD monitors

Merivaara offers a wide range of medically certified monitors suitable for visualizing the content during medical procedures. We offer also a range of compatible accessories and even customized mounting solutions to meet your demands.



### Available models

#### **SUPPLIER CODE**

FS-L2403D  
FS-L2702D

#### **PRODUCT CODE**

A45000185  
A45000472

## **Contact us**

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)