## **SOLO**

# SOLO suspension solution for surgical lights and monitors

Modern operating rooms today are equipped with more and more equipment designed to provide surgical staff both the clinical images as well as critical video sources during the procedure. For the staff, functionality and ergonomics of suspension solutions that carry these vital monitors, is of very high importance.

Suspension solutions of Q-Flow product family has gone through a complete update, now being able to offer many new solutions that will enhance the ergonomics of any operating room.

Read how Merivaara's R&D was challenged to think outside the box and what then happened.

# **EXAMPLES OF AVAILABLE ARM OPTIONS:**

- SOLO Mobile
- SOLO Wall
- **SOLO Ceiling** 50 mm, 65 mm or 110 mm ceiling tube







Merivaara Corp. 24.11.2025



### **Technical Specifications**

Product Code	520210

#### **Contact us**

Merivaara Corp. Tarmontie 2-4 15860 HOLLOLA FINLAND

+358 3 3394 611

merivaara@merivaara.com

Merivaara Corp. 24.11.2025