


## Back and kidney elevator 19280

Back and kidney elevator 19280 to be placed on the table top or on the mattress for kidney surgical procedures. The angle is adjustable with a crank mechanism, 245 degrees. A separate VEF mattress is included 500 x 1080 mm.



### Technical Specifications

Length	555 mm
Height	540 mm
Height adjustment	60 - 190 mm
Max. patient weight	180 kg
Safe working load (SWL)	75 kg
Weight	8.8 kg
Product Code	100019280

## **Contact us**

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)