


## Oxygen bottle holder 18189

For holding 5-10 L oxygen bottle on the Rapido operating tables's base cover.

- For bottle diameter ca. 140 mm
- Max bottle length with valves 550 mm
- To be fitted by the customer, mounting instructions included
- Nylon fabric straps with zinc plated locking buckles
- Suitable with Rapido



### Technical Specifications

Вес	2.25 kg
код продукта	100018189

### Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)