

Arm and hand table rectangular AB5010

- light and user-friendly arm and hand surgery table with two Easy clamps CL2010
- table top enables fluoroscopy, X-ray imaging area in the middle of table 350 mm
- height adjustable
- rod: 135 mm
- includes mattress AB5101 Supporting leg





TECHNICAL SPECIFICATIONS

| | |
|--------------|---------------------------------------|
| | |
| Length | 650 mm |
| Width | 400 mm |
| Product code | AB5010 (EU), AB5011 (US), AB5012 (UK) |

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com