

Arm rest 125

Padded arm board 125 includes a fixing clamp and straps 110. Adjustable height, gas spring assisted angle adjustment. Rotates 180 degrees. Suitable for EUR, US and UK rails.



Technical Specifications

Width	640 mm
Height	205 mm
Thickness	60 mm
Max. patient weight	325 kg
Safe working load (SWL)	19.5 kg
Trendelenburg	-36 °
Reverse Trendelenburg	+26 °
Weight	6.2 kg

Product Code

110125

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com