Arm rests 151

Arm rest pair 151 for opthalmology is designed to support surgeon's arms during operation. Arm board turns 360 degrees. Height adjustment with joint screw on the side of the aluminium clamp. Fits all Merivaara operating tables. Rod diameter 16 mm. Including clamps for fixing. Suitable for EUR, US and UK rails.





Technical Specifications

Width	300 mm
Length	85 mm
Height adjustment	120 - 320 mm
Thickness	40 mm q
Safe working load (SWL)	10 kg
Weight	2.4 kg
Product Code	100000151

Merivaara Corp. 26.10.2025



Contact us

Merivaara Corp. Tarmontie 2-4 15860 HOLLOLA FINLAND

+358 3 3394 611

merivaara@merivaara.com

Merivaara Corp. 26.10.2025