

Armboard AB2110

- one-handed adjustment armboard with Easy clamp
- wide range of height and angle adjustment with ball joint
- antimicrobial adjustment handle
- rod: 150 mm
- includes VEF mattress



TECHNICAL SPECIFICATIONS

Length	620 mm
Width	170 mm
Max. patient weight	250 kg
Safe working load (SWL)	21 kg
Longitudinal tilt (37,5°/direction)	75°
Lateral tilt (37,5°/direction)	75°
Sideways movement (rotates around a ball joint)	360°
Weight	4,4 kg
Product code	AB2110 (EU), AB2111 (US), AB2112 (UK)

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com