

## Armboard long rod AB2130

- long rod model
- one-handed adjustment armboard with Easy clamp
- wide range of height and angle adjustment with ball joint
- antimicrobial adjustment handle
- rod: 280 mm
- includes VEF mattress



### TECHNICAL SPECIFICATIONS

Length	620 mm
Width	170 mm
Max. patient weight	250 kg
Safe working load (SWL)	21 kg
Longitudinal tilt (37,5°/direction)	75°
Lateral tilt (37,5°/direction)	75°
Sideways movement (rotates around a ball joint)	360°
Weight	4,6 kg
Product code	AB2130 (EU), AB2131 (US), AB2132 (UK)

## **Contact us**

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)