

## HPL-arm rest 60510

HPL-arm rest 60510 is specially made to be used with the carbon fibre back rest. The patient's counterweight holds the arm board in place. No metal parts. Total length 900 mm, arm rest size 160 x 680 mm. Antistatic VEF mattress of 80 mm thickness included.



### Technical Specifications

Product code	100060510

## **Contact us**

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)