


# Multi-task armboard MV-11802

Arm rest provides arm support for prone, neuro, lateral and chair positions. A wide range of position flexibility, turns 360 degrees and is height adjustable. The premium durasoft pad ensures patient comfort. Universal clamp 10300 is needed for fixing.

## Order no.

110802	Multi-task armboard MV-11802 with
100 011 802	Lateral arm board MV-11802
100 010 300	Universal clamp 10300



## Technical Specifications

Width	510 mm
Length	150 mm
Height adjustment	170 - 480 mm
Thickness	30 mm
Max. patient weight	225 kg
Safe working load (SWL)	16 kg
Weight	2.2 kg
Lateral tilt	45 - 45 °

Product Code

110802, 100011802

## Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)