

Easy clamp Ø18 mm rotation

- light and user-friendly aluminium rotation lever clamp
- rotation feature allows versatile functionality and adjustability
- antimicrobial tightening levers
- for attaching rods with a diameter of 16-18 mm



TECHNICAL SPECIFICATIONS

Product code	CL2110 (EU), CL2111 (US), CL2112 (UK)

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com