

Extension rail 19128

Extension rail 19128 for extending the length of the table side rails. Fits to both sides of the table. Stainless steel.

Order no.

110128 Rail extension, incl. clamp for US or UK rails

100 019 128U (For US rail)



Technical Specifications

| | |
|-------------------------|--------------------|
| | |
| Length | 65 mm |
| Height | 550 mm |
| Max. patient weight | 225 kg |
| Safe working load (SWL) | 58 kg |
| Weight | 2.8 kg |
| Product Code | 110128, 100019128U |

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com