

## Transfer mattress 19400

Used for transferring a patient between different platforms. Easy to use, easy to wipe clean, waterproof material. Colour dark blue. CE marked. Not intended for lifting the patients.



### Technical Specifications

Width	550 mm
Length	1800 mm
Thickness	30 mm
Weight	1.9 kg
Product Code	100019400

### Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)

